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1. Allen KD, Larzelere RE, Hutfless S, Beiraghi S. Comparison of a computerized anesthesia device with a traditional syringe in preschool children. *Pediatr Dent.* 2002;24:315-320.
2. Ram D, Peretz B. Assessing the pain reaction of children receiving periodontal ligament anesthesia using a computerized device (Wand). *J Clin Pediatr Dent.* 2003;27:247-250.
3. Ashkenazi M, Blumer S, Ilana E. Effectiveness of computerized delivery of intrasulcular anesthetic in primary molars. *JADA.* 2005;136:1418-1425.
4. Versloot J, Veerkamp JSJ, Hoogstraten J. Computerized anesthesia delivery system vs. traditional syringe: comparing pain and pain-related behavior in children. *Eur J Oral Sci.* 2005;113:488-493.
5. Öztas N, Ulusu T, Bodur H, Dogan C. The Wand in pulp therapy : An alternative to inferior alveolar nerve block. *Quint. International.* 2005;36:(7)559-564.
6. Nicholson JW, Berry TG, Summitt JB, Yuan CH, Witten TM. Pain perception and utility: A comparison of the syringe and computerized local injection techniques. *Gen Dent.* 2001;167-172.

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⁵ The AMSA injection: A new concept for local anesthesia of maxillary teeth using a computer-controlled injection system; Friedman, Mark J., DDS; Hochman, Mark N., DDS; Quintessence Int. 1998: 29:297-303.

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1. Hochman MN, Chiarello D, Hochman CB, Lopatkin R, Pergola S. Computerized Local Anesthesia Delivery vs. Traditional Syringe Technique. NY State Dent J. 1997;63:24-29.
2. Friedman MJ, Hochman MN. 21st Century Computerized Injection for Local Pain Control. Compend Contin Educ Dent. 1997;18:995-1003.
3. Krochak M, Friedman N. Using a precision-metered injection system to minimize dental injection anxiety. 1998;19(2):137-148.
4. Friedman MJ, Hochman MN. The AMSA injection: A new concept for local anesthesia of maxillary teeth using a computer-controlled injection system. Quintessence Int. 1998; 29:297-303.
5. Farah JW. Editors Choice-The Wand. The Dental Advisor. 1998;15:1.
6. CRA. Local Anesthesia, Automated Delivery. Clinical Research Associates Newsltr. 1998;22:1-2.
7. Friedman MJ, Hochman MN. P-ASA Block Injection: A new palatal technique to anesthetize maxillary anterior teeth. J Esthet Dent. 1999;11(2):63-71.
8. Gardner M. The AMSA Block: It will leave your patients smiling. Oral Health. 1999;July:43.
9. Leiberman, William H. Clinical Session:The Wand. Pediatric Dent. 1999;21:2.
10. Levato C. Giving the Wand a shot. Dent Pract Fin. 1998;July:53.
11. Kehoe B. In search of the painless injection. Dent Pract Fin. 1998;July:53.
12. Lackey A. Technology:An advancement in the delivery of local anesthesia. Pract Perio Aesthet Dent. 1998;10:1191-1193.
13. Asarch T, Allen K, Petersen B, Beiraghi S. Efficacy of a computerized local anesthesia device in pediatric dentistry. 1999;21:421-424.
14. Oldak S, Jackson LA. As we see it: The Wand. J Southeast Soc Pediat Dent.

1999;5:38.

15. Kronish E. Creating a less painful image of dentist. AGD Impact. 1999;April
16. Hochman MN, Friedman MJ. In vitro study of needle deflection: A linear insertion technique versus a bi-directional rotation insertion technique. Quintessence Int. 2000;31:737-743.
17. Gibson RS, Allen K, Hutfless S, Beiraghi S. The Wand vs. traditional injection: A comparison of pain related behaviors. Pediatric Dent. 2000;22:458-462.
18. Froum SJ, Tarnow D, Caiazzo A, Hochman MN. Histologic response to intraligament injections using a computerized local anesthetic delivery system. A pilot study in Mini-Swine. J Periodontol. 2000;71:1453-59.
19. Lipton L. Using Computer-controlled technology to alleviate stress & reduce discomfort during local anesthetic delivery in a pediatric practice. J Southeastern Soc Pediatric Dent. 2000;6:22-32
20. Goodell GG, Gallagher FJ, Nicol BK. Comparison of a controlled injection pressure system with a conventional technique. Oral Surg Oral med Oral Pathol Oral Radiol Endod. 2000;90:88-94.
21. Friedman MJ, Donaldson D, Malamed SF, Yagiela JA. Technology Forum: New Advances in Local Anesthesia. Compend Contin Educ Dent. 2000;21:432-440.
22. Grace EG, Barnes DM, Macek MD. Patient and dentist satisfaction with a computerized local anesthetic injection system. Compend Contin Educ Dent. 2000;21:746-752.
23. Aboushala A, Kugel G, Efthimiadis N, Korchak M. Efficacy of a computer-controlled injection system of local anesthesia in vivo. IADR Abstract. 2000;Abst#2775.
24. Cheng H, Pong PY, Chang WJ, Lee SY. Using a computer-controlled injection system to minimize dental injection pain. IADR Abstract. 2000;Abst#2777.
25. Loomer PM, Perry DA. Efficacy of computer-controlled local anesthesia during scaling and root planing. IADR Abstract. 2000;Abst#590.
26. Koili K, Boyles J, Gavlak J, Weaden S, Crout R. Comparing the efficacy of the Wand and traditional buccal infiltrations. IADR Abstract. 2000;Abst#2772.
27. Nicholson JW, Berry TG, Summitt JB, Yuan CH, Witten TM. Pain perception and utility: A comparison of the syringe and computerized local injection techniques. Gen Dent. 2001;167-172.

28. Jackman DS, Hertz MB. Techniques of Drug Administration. Oral Maxillo Surg Clinics North Amer. 2001;13:199-213.
29. Friedman MJ, Hochman MN. Using AMSA and P-ASA nerve blocks for esthetic restorative dentistry. Gen Dent. 2001;49(5):506-511.
30. Hochman MN, Friedman MJ. An in vitro study of needle force penetration comparing a standard linear insertion to the new bidirectional rotation insertion technique. Quintessence Int. 2001;32:789-796.
31. Fukayama H. New Trends in Local Anesthesia. Hyogo Dental Assoc J. 2001;Jan;593-602.
32. Tan PY, Vukasin P, Chin ID, Ciona CJ, Ortega AE, Anthone GJ, Corman ML, Beart RW. The Wand local anesthetic delivery system. Diseases Colon Rectum. 2001;44:686-689.
33. Landsman A, DeFronzo D, Hedman J, McDonald J. A new system for decreasing the level of injection pain associated with local anesthesia of a toe. Am Acad Podiat Med. 2001;Abstract.
34. Barusco MN, Leavitt ML. The use of computerized anesthesia injection system to minimize pain during hair transplant surgery. Hair Transplant Forum Inter. 2001;11:107-108.
35. Isen D. A review of computer controlled injection devices. Oral Health. 2001 July:31-34.
36. Kudo M, Ohke H, Katagiri K, Sato Y, Kawai T, Kato M, Kokubu M, Shinya N. The shape of local anesthetic injection syringes with less discomfort and anxiety- Evaluation of discomfort and anxiety caused by various types of local anesthetic injection syringes in high level trait-anxiety people. J Japan Dent Soc Anesthesiol. 2001;29:173-178.
37. Rosenberg E. A computer-controlled anesthetic delivery system in a periodontal practice: Patient satisfaction and acceptance. J Esthet Restor Dent. 2001;13:25-32.
38. Allen KD, Kotil D, Larzelere RE, Hutfless S, Beiraghi S. Comparison of a computerized anesthesia device with a traditional syringe in preschool children. Pediatr Dent. 2002 Jul-Aug;24(4):315-20.
39. True RH, Elliot RM. Microprocessor-controlled local anesthesia versus the conventional syringe technique in hair transplantation. Dermatol Surg. 2002;28:64-69.

40. Swanepoel PF, Heystek P. Computer assisted local anesthetic application for nasal surgery. 8th AAFPRS Inter. Sympos. 2002;Abstract.
41. JADA. Dental Product Spotlight:Local anesthetic delivery. JADA 2002;133(JADA's 1st product review);106.
42. Blanton PL, Jeske AH. Dental Local Anesthetics: Alternative Delivery Methods. 2003;134:228-234.
43. Perry DA, Loomer PM. Maximizing Pain Control. The AMSA Injection can provide anesthesia with few injections and less pain. Dimensions of Dental Hygiene 2003;April/May:28-33.
44. Ram D, Peretz B .J Clin Pediatr Dent. 2003 Spring;27(3):247-50.Assessing the pain reaction of children receiving periodontal ligament anesthesia using a computerized device (Wand).
45. Fukayama H, Yoshikawa F, Kohase H, Umino M, Suzuki N. Efficacy of AMSA anesthesia using a new injection system, the Wand. Quintessence International, 2003;34:537-541.
46. Peter M. Loomer & Dorothy A. Perry, Comparison of Computer-Controlled Delivery to Syringe Delivery of Local Anesthetic During Therapeutic Scaling and Root Planing. JADA 2004;135:358-365.
47. Kasaj A, Berakdar M, Nicolaescu A, Willershausen, Sculean A. Evaluation of a new anesthesia technique for nonsurgical periodontal therapy. Johannes Gutenberg-University, Mainz, Germany. IADR/AADR/CADR 82nd General Session (March 10-13, 2004) Abstract #222 - Anesthesiology Research 2
48. Schwartz-Arad D, Dolev E, Williams W. Maxillary nerve block – A new approach using a computer-controlled anesthetic delivery system for maxillary sinus elevation procedure. A prospective study. Quintessence International, 2004;35:477-480.
49. CRA Newsletter. Products reported most by CRA evaluators. Products CRA evaluators “Can’t Live Without”. July 2004.;28(7):2-4.
50. Palm AM, Kirkegaard U, Poulsen S. The Wand versus Traditional Injection for Mandibular Nerve Block in Children and Adolescents: Perceived Pain and Time of Onset. Pediatric Dentistry, 2004;26:481-484.

51. Shepherd PA, Eleaszer PD, Clark SJ, Scheetz JP. Measurement of Intraosseous Pressures Generated by the Wand[®], High-Pressure Periodontal Ligament Syringe, and the Stabident System. *J. Endodontics*, 2001;27(6):381-384.
52. Ashkenazi M, Blumer S, Eli I. Effective of Computerized Delivery of Intrasulcular Anesthetic in Primary Molars. *JADA*, 2005;136:1418-1425.
53. Ghelber O, Gebhard R, Adebayo G, Szmuk P, Hagberg C, Ianucci D.:Utilization of the CompuFlo™ in determining the pressure of the epidural space: a pilot study. *Anesth Analg* 2005;100:S-189.
54. Ghelber O, Gebhard R, Szmuk P, Hagberg C, Ianucci D.: Identification of the epidural space-a pilot study of a new technique. *Anesth Analg* 2005;100:S-255
55. Gebhard R, Ghelber O, Szmuk P, Pivalizza E, Walters D: Pressure Monitoring During Injection of Local Anesthetics for Nerve Blocks Utilizing the CompuFlo[®] Injection Pump. *Anesth Analg* 2005
56. Kudo M. Initial Injection Pressure for Dental Local Anesthesia: Effects on Pain and Anxiety. *Anesthesia Progress*, 2005; 52:95-101.
57. Versloot J, Veerkamp JSJ, Hoogstraten J. Computerized anesthesia delivery system vs. traditional syringe: comparing pain and pain-related behavior in children. *Eur J Oral Sci.* 2005;**113**:488-493.
58. Öztas N. Ulusu T. Bodur H. Dogan C. The Wand in pulp therapy : An alternative to inferior alveolar nerve block. *Quint. International.* 2005;**36**:(7)559-564.
59. Ram D. Kasssire J. Assessment of a palatal approach-anterior superior alveolar (P-ASA) nerve block with the Wand[®] in paediatric dental patients. *J Clin Pediatr Dent.* 2006;**16**:348-3551.
60. Jalevik B, Klingberg G, and G. KLINGBERG, Sensation of Pain when using Computerized Injection Technique, the Wand™. *IADR Pan European Federation Sept. 13-16, 2006.*
61. Hochman MN, Friedman MF, Williams WP, Hochman CB. Interstitial Pressure Associated with Dental Injections: A Clinical Study. *Quintessence International*, 2006;37:469-476.

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1. Michaelian MJ, Agha-razi F, Hutter J. Anesthetic efficacy of the periodontal ligament injection using the Wand vs. the intra-osseous injection using stabident. (Unpublished manuscript, BU Dental)

2. Franco L, Naseri L, Hochman MN, Camarda AJ. A New Multi- Cartridge Injection Technique for Achieving Safe and Effective Dental Local Anesthesia. Submitted for publication, Oct. 2003.

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Barnard D. Hazards of Local Anesthesia Injections. ISBN: 0-620-26308-3. Pretoria, South Africa. (pg2) 1998.

Murphy D. Ergonomics and the Dental Care Worker. ISBN: 0-87553-0233-0. Washington DC, American Public Health Association. (pg 181) 1998.

Wilkins E. Clinical Practice of the Dental Hygienist 8th Ed. ISBN: 0-683-30362-7. Philadelphia, Pennsylvania. (pg 503) 1999.

Dionne R, Phero J, Becker D. Management of Pain and Anxiety in the Dental Office. ISBN: 0-7216-7278-7. Philadelphia, Pennsylvania. (pg 204-06) 2002.

Malamed S. Handbook of Local Anesthesia 5th Ed. ISBN:0-323-02449-1. Elsevier/Mosby, St. Louis, Missouri. 2004.

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©2007-2023 Milestone Scientific Inc.