

MILESTONE
SCIENTIFIC®

COMPUFLO®
INSTRUMENTS

The Modern Epidural Solution





WE HAVE THE SOLUTION

CompuFlo® Instrument

A unique opportunity to add **efficiency, safety** and **satisfaction** with the CompuFlo Epidural Instrument – a modern solution to objectively identify the epidural space with a **99% success rate** on the first attempt.

WHAT IS THE SYSTEM?

The CompuFlo Epidural Instrument allows for a more precise, reliable and effective application of epidural anesthesia through real time needle location. The instrument offers objective detection of pressure changes to confirm an epidural space within minutes.

HOW DOES THE SYSTEM WORK?

Our patented DPS Dynamic Pressure Sensing Technology® monitors the exact pressure, and optimal needle position during the administration of anesthetic. Visual and audible feedback from the device, aids in the identification of the correct needle and catheter placement for labor and delivery epidural procedures.

DID YOU KNOW?

- ✓ 82% of anesthesiologists believe there are hospital mandated time pressures that cause unnecessary risks
- ✓ 58% are concerned about dural punctures during labor and delivery epidural procedures
- ✓ 96% believe that they will need to take on more of a leadership role in making decisions that will impact patient care

*Research was conducted by Schlesinger Group, September 2020

“The CompuFlo Epidural Instrument improves the block. It does help you to make a decision and decrease the complications; reducing the length of stay and improve patient satisfaction.”

– Dr. Ayman A Alian, MBChB,
Associate Professor of Anesthesiology Division Chief,
Obstetric & Gynecological Anesthesiology – Yale School of Medicine


WHAT ARE THE BENEFITS

CompuFlo Epidural Instrument

- ✓ Increase patient safety by avoiding dural punctures
- ✓ More than **99%** effective on the first attempt
- ✓ Shortening the training time of an anesthesiologist resident
- ✓ Avoid the cost of litigation
- ✓ Save an average of **\$504** per hospital stay with CompuFlo Epidural Instrument compared to the traditional LOR technique

EFFICIENCY, SAFETY AND COMFORT

- ✓ 25% reduction in procedure time for labor and delivery
- ✓ Detects subtle pressure changes four times a second
- ✓ Immediate catheter verification with CathCheck™ Verification System within one to two minutes
- ✓ Precisely controls the flow rate of fluid with real-time feedback



"In the long run, the CompuFlo Epidural Instrument will help a lot of patients, helping them be more comfortable and prevent the side effects of epidurals."

- Dr. Abraham Scheer

“Our vision is to take drug
delivery from the 1900s
to the 21st century –
now it’s time to take it further.”

How Do I Get Started?

1-800-862-1125

COMPUFLO[®]
INSTRUMENTS

VISIT US ONLINE

milestonescientific.com

Milestone Scientific Inc. • 425 Eagle Rock Avenue, Suite 403 Roseland, NJ 07068

ML-0270 Rev. A